

IR Spectroscopy Analyzer

Group 12

Asia 黃夏彤 Charoline 黃菊花 Teresa 陳翠珍 Thalia 郭伶美

Basic concepts

What is IR spectroscopy?


A method used in analytical chemistry to predict a chemical structure using infrared electromagnetic waves

Usually, IR is used to confirm the existence of a functional group for chemical synthesis!

How it works?


Different functional groups and bonds affect the waves' absorbance at various wavelengths. This data is then plotted and used to analyze the compound's structure!

Functional group	Mode of vibration	Wavenumber (cm^{-1})	Intensity
C-H	Alkanes (stretching)	3000-2850	Strong
	-CH ₃ (bending)	1450 and 1375	Medium
	-CH ₂ - (bending)	1465	Medium
Alkenes	(stretching)	3100-3000	Medium
	(out of plane bending)	1000-650	Strong
	Aromatics (stretching)	3150-3050	Strong
Alkyne	(out of plane bending)	900-690	Strong
	(stretching)	ca. 3300	Strong
	Aldehyde	2900-2800 2800-2700	Weak Weak
C=C	Alkane	not interpretative useful	
	Alkene	1680-1600	Medium to Weak
	Aromatic	1600 and 1475	Medium to Weak
C≡C	Alkyne	2250-2100	Medium to Weak



Why?

- Manually pick out which peak to use -> bias
- Sometimes peaks are too close -> hard to notice
- Checking the table now and then is time-consuming





Preprocessing Data

The original data from IR spectroscopy is typically presented in terms of transmittance, which can contain significant biases. Therefore, the raw data is processed using the equation

$$A = -\log_{10}(T)$$

Purposes:


- 1. Transforms the relationship:** It converts transmittance (which can decrease in a valley-like shape) into absorbance
- 2. Enhances clarity:** The logarithmic function magnifies differences in low transmittance values, making small features or variations in the data more visible and distinct




Preprocessing Data

$$A = -\log_{10}(T)$$

Relationship of Absorbance and Transmittance





CODE

Code Overview

File Conversion

JD^X FILE

X1 Y1 Y2 Y3 Y4 Y5
X6 Y6 Y7 Y8 Y9 Y10

.....



CSV FILE

X1 Y1
X2 Y2
X3 Y3

...

Peak Detection

Use Mean and STD
to find threshold



Use `find_peaks`
from `scipy`



Handle overlapping
peaks

Categorization


Create a dictionary
for functional
groups & categorize
the peaks



Categorize the
intensity

Printed Results

Create plots with
matplotlib, tables
with tabulate, and
interactive display
with ipywidgets



Data Source : NIST Chemistry WebBook, SRD 69

File Conversion

JD^X FILE

X1 Y1 Y2 Y3 Y4 Y5

X6 Y6 Y7 Y8 Y9 Y10




CSV FILE

X1 Y1

X2 Y2

X3 Y3

...



```
##MAXY=1.05
##MINY=0.00216665
##NPOINTS=3576
##XYDATA=(X++(Y..Y))
587.635000 0.7290 0.7270 0.7240 0.7200 0.7110
592.271792 0.7060 0.6950 0.6840 0.6681 0.6581
596.908584 0.6460 0.6450 0.6460 0.6470 0.6530
601.545377 0.6570 0.6600 0.6620 0.6630 0.6680
606.182169 0.6740 0.6760 0.6840 0.6910 0.6940
610.818961 0.6970 0.7050 0.7060 0.7070 0.7080
615.455753 0.7120 0.7190 0.7220 0.7240 0.7270
620.092545 0.7310 0.7380 0.7410 0.7420 0.7450
624.729337 0.7490 0.7520 0.7580 0.7580 0.7600
```

```

import csv
import math
def read_jdx_file(file_path):

#Reads a JCAMP-DX file, extracts the X and Y values, and returns a list of (x, y) pairs.
    with open(file_path, 'r') as file:
        lines = file.readlines()

    x_values = []
    y_values = []
    superX_values=[]
    superY_values=[]

    # Skip all lines that start with '#' and begin processing the data
    data_started = False
    for line in lines:
        line = line.strip()

        # Skip lines that did not start with a number
        if not line or not line[0].isdigit():
            continue

        # If we encounter a non-# line, data starts
        data_started = True

        # Process data lines only
        if data_started and line:
            # Split the line into individual values
            values = list(map(float, line.split()))

            superX_values.append(values[0]) # First value is X
            superY_values.append(values[1:]) # Remaining values are Y
#result -> superX_values=[x1, x6, ...] superY_values=[(y1,y2,y3,y4,y5), (y6, y7, y8, y9, y10), ...]

    for x in superX_values:
        x_values.append(x)
        x_values.append(x+1)
        x_values.append(x+2)
        x_values.append(x+3)
        x_values.append(x+4)
#result -> x_values=[x1, x2, x3,...]

    for ys in superY_values:
        for y in ys:
            y_values.append(y)
#result -> y_values=[y1, y2, y3,...]

    if not x_values or not y_values:
        raise ValueError("No data found in the JCAMP-DX file.")

```

Read JDX file

Create empty list for X and Y Values

Skips line with no data

Added the missing increment

X1 X5 X10 [superX_values] -> X1 X2 X3... [x_values]

Convert

(Y1,Y2...,Y5), (Y6,...,Y10)... [superY_values] -> Y1, Y2, Y3... [y_values]

File Conversion

JDX FILE

X1 Y1 Y2 Y3 Y4 Y5

X6 Y6 Y7 Y8 Y9 Y10

.....



CSV FILE

X1 Y1

X2 Y2

X3 Y3

...

```

def generate_xy_pairs(x_values, y_values):
    """
    Converts the list of X values and corresponding list of Y values into (x, y) pairs.
    """
    xy_pairs = []
    counter=0
    while counter < len(y_values):
        y=y_values[counter]
        x=x_values[counter]
        if y==0:
            xy_pairs.append((x, y))
        else:
            y=-math.log(y,10)
            xy_pairs.append((x, y))
        counter+=1
    return xy_pairs
#result: [(x1, y1), (x2, y2), (x3, y3), ...]

```

```

def save_xy_to_csv(output_path, xy_pairs):

#Saves the (x, y) pairs to a CSV file.
    with open(output_path, mode='w', newline='') as file:
        writer = csv.writer(file)

        # Write the header
        writer.writerow(["Wavenumber", "Absorbance"])

        # Write the (x, y) pairs
        for x, y in xy_pairs:
            writer.writerow([x, y])

```

**Output: A CSV file with two columns:
Wavenumber and Absorbance**

**Generate (x,y) pairs from the extracted
x_values and y_values**

$$A = -\log_{10}(T)$$

Save the (x, y) pairs to a CSV file

```

# Main function
if __name__ == "__main__":
    input_file = "formaldehyde.jdx" # Replace with the path to your JCAMP-DX file
    output_file = "formaldehyde.csv" # Output file for (x, y) pairs in CSV format

    try:
        # Read the .jdx file, skipping the first 36 lines
        x_values, y_values = read_jdx_file(input_file)

        # Generate the (x, y) pairs
        xy_pairs = generate_xy_pairs(x_values, y_values)

        # Save the results to a .csv file
        save_xy_to_csv(output_file, xy_pairs)

        print(f"Converted (x, y) data saved to {output_file}")
    except Exception as e:
        print(f"Error: {e}")

```

Converted (x, y) data saved to formaldehyde.csv

```

import numpy as np
from ipywidgets import interact, widgets
import pandas as pd
from scipy.signal import find_peaks
import matplotlib.pyplot as plt
from tabulate import tabulate
from IPython.display import display

# Load data
data = pd.read_csv('formaldehyde.csv')

# Check for peaks in 'Absorbance'
# Calculate the mean and standard deviation for threshold
mean_height = data['Absorbance'].mean()
std_height = data['Absorbance'].std()
threshold = mean_height + std_height

# Find peaks above the threshold
peak_indices, properties = find_peaks(data['Absorbance'], height=threshold)

# Loop through all the peaks and check for nearby peaks
unique_peak_indices = []
max_distance = 15 # The distance for overlapping peaks

for i in peak_indices:
    is_unique = True
    for j in unique_peak_indices:
        if abs(data['Wavenumber'][i] - data['Wavenumber'][j]) < max_distance:
            if data['Absorbance'][i] > data['Absorbance'][j]:
                unique_peak_indices.remove(j)
                unique_peak_indices.append(i)
            is_unique = False
            break
    if is_unique:
        unique_peak_indices.append(i)

```

Peak Detection

Use Mean and STD
to find threshold



Use `find_peaks`
from `scipy`



Handle overlapping
peaks

```

# Dictionary for functional groups
functional_groups = {
    "C-H Alkanes": (3000, 2850),
    "C-H Alkanes -CH3": [(1450, 1450), (1375, 1375)],
    "C-H Alkanes -CH2-": (1465, 1465),
    "C-H Alkenes": (3100, 3000),
    "C-H Aromatics": [(3150, 3050), (900, 690)],
    "C-H Alkyne": (3300, 3300),
    "C-H Aldehyde": [(2900, 2800), (2800, 2700)],
    "C=C Alkene": (1680, 1600),
    "C=C Aromatic": [(1600, 1600), (1475, 1475)],
    "C≡C Alkyne": (2250, 2100),
    "C=O Aldehyde": (1740, 1720),
    "C=O Ketone": (1725, 1705),
    "C=O Carboxylic acid": [(1725, 1700), (3400, 2400)],
    "C=O Ester": (1750, 1730),
    "C=O Amide": (1670, 1640),
    "C=O Anhydride": [(1810, 1810), (1760, 1760)],
    "C=O Acid chloride": (1800, 1800),
    "C-O Alcohols, Ethers, Esters": (1300, 1000),
    "O-H Alcohols, Phenols (Free)": (3650, 3600),
    "O-H Alcohols, Phenols (H-Bonded)": (3500, 3200),
    "N-H Primary and Secondary Amines": [(3500, 3100), (1640, 1550)],
    "C-N Amines": (1350, 1000),
    "C=N Imines and Oximes": (1690, 1640),
    "C≡N Nitriles": (2260, 2240),
    "X=C=Y Allenes, Ketenes, Isocyanates": (2270, 1950),
    "N=O Nitro (R-N02)": [(1550, 1550), (1350, 1350)],
    "S-H Mercaptans": (2550, 2550),
    "S=O Sulfoxides": (1050, 1050),
    "C-X Fluoride": (1400, 1000),
    "C-X Chloride": (800, 600),
    "C-X Bromide and Iodide": (667, 0),
}

```

Categorization

Create a dictionary for functional groups



Create and design the scatter plot


Printed Results

```

# Create label: 1 if there's a peak, otherwise 0
data['label'] = 0
data.loc[unique_peak_indices, 'label'] = 1

# Visualize data to verify peaks
plt.figure(figsize=(10, 6))
plt.plot(data['Wavenumber'], data['Absorbance'], label='Absorbance')
plt.scatter(data.loc[unique_peak_indices, 'Wavenumber'], data.loc[unique_peak_indices, 'Absorbance'], color='red', label='Peaks')
plt.xlabel('Wavenumber')
plt.ylabel('Absorbance')
plt.title('Absorbance vs Wavenumber with Peak Detection')
plt.legend()
plt.show()
print("\n")

```



```

def combined_analysis(start_wavenumber=None, end_wavenumber=None):
    print("\n")

    # Filter data based on the user-defined range
    if start_wavenumber and end_wavenumber:
        filtered_data = data[(data['Wavenumber'] >= start_wavenumber) & (data['Wavenumber'] <= end_wavenumber)]
        unique_peak_indices_filtered = [i for i in unique_peak_indices if start_wavenumber <= data['Wavenumber'][i] <= end_wavenumber]
    else:
        filtered_data = data
        unique_peak_indices_filtered = unique_peak_indices

    if len(unique_peak_indices_filtered) == 0:
        print("No peaks detected in the selected range. Please select another range.")
        return # Exit the function if no peaks are found

    peak_values = filtered_data.loc[unique_peak_indices_filtered, ['Wavenumber', 'Absorbance']]
    peak_values['Wavenumber'] = peak_values['Wavenumber'].astype(float)
    peak_values['Absorbance'] = peak_values['Absorbance'].astype(float)

    # Identify the highest peak
    highest_peak_index = peak_values['Absorbance'].idxmax()
    highest_peak_row = peak_values.loc[highest_peak_index]

    # Prepare rows for the combined table
    table_data = []

    # Header row for the table
    table_data.append(["Wavenumber", "Absorbance", "Functional Groups", "Peak Intensity"])

    # Define thresholds for categorization of high, medium, and low
    peak_heights = peak_values['Absorbance'].values
    high_threshold = np.percentile(peak_heights, 80) # Top 20%
    low_threshold = np.percentile(peak_heights, 20) # Bottom 20%

```

Make sure that the input from the user does not result in an error

Change the type of the data

Identify the peaks

Create the table for the printed data

Define the threshold for the intensity of the peaks

Preventing out of range error

```
# Filter data based on the user-defined range
if start_wavenumber and end_wavenumber:
    filtered_data = data[(data['Wavenumber'] >= start_wavenumber) & (data['Wavenumber'] <= end_wavenumber)]
    unique_peak_indices_filtered = [i for i in unique_peak_indices if start_wavenumber <= data['Wavenumber'][i] <= end_wavenumber]
else:
    filtered_data = data
    unique_peak_indices_filtered = unique_peak_indices

if len(unique_peak_indices_filtered) == 0:
    print("No peaks detected in the selected range. Please select another range.")
    return # Exit the function if no peaks are found
```

Default range: Start Wavenumber = 587.63, End Wavenumber = 3902.94
If you wish to analyze a specific range, you can enter the range below:

Start Wave...

End Waven...

No peaks detected in the selected range. Please select another range.

Creating Printed Tables

```
# Prepare rows for the combined table
table_data = []

# Header row for the table
table_data.append(["Wavenumber", "Absorbance", "Functional Groups", "Peak Intensity"])
```

Wavenumber	Functional Groups	Absorbance

```

# Loop through each peak to add functional groups and intensity categorization
for index, row in peak_values.iterrows():
    wavenumber = row['Wavenumber']
    absorbance = row['Absorbance']
    matched_groups = []

    # Find matching functional groups
    for group, value in functional_groups.items():
        if isinstance(value[0], int): # Single range
            high, low = value
            if low <= wavenumber <= high:
                matched_groups.append(group)
        else: # Multiple ranges
            for high, low in value:
                if low <= wavenumber <= high:
                    matched_groups.append(group)

    # Categorize peak intensity
    category = "High" if absorbance >= high_threshold else "Low" if absorbance <= low_threshold else "Medium"

    # Style the highest peak row
    if index == highest_peak_index:
        wavenumber = f"\033[1m\033[31m{wavenumber:.2f}\033[0m" # Bold + Red
        category = f"\033[1m\033[31m{category}\033[0m" # Bold + Red

    # Add row to the table
    functional_groups_str = ', '.join(matched_groups) if matched_groups else "No matching functional group"
    table_data.append([wavenumber, functional_groups_str, category])

```

Create a for loop for repetition

Check if the peak is in the range of the functional group

Categorize the intensity of each peak



Label the highest peak in red

Print out the result

Categorizing the Intensity

```
# Define thresholds for categorization of high, medium, and low
peak_heights = peak_values['Absorbance'].values
high_threshold = np.percentile(peak_heights, 80) # Top 20%
low_threshold = np.percentile(peak_heights, 20) # Bottom 20%

# Loop through each peak to add functional groups and intensity categorization
for index, row in peak_values.iterrows():
    # Categorize peak intensity
    category = "High" if absorbance >= high_threshold else "Low" if absorbance <= low_threshold else "Medium"
```



Matching the Peaks to the Functional Groups

```

for index, row in peak_values.iterrows():
    wavenumber = row['Wavenumber']
    absorbance = row['Absorbance']
    matched_groups = []

    # Find matching functional groups
    for group, value in functional_groups.items():
        if isinstance(value[0], int): # Single range
            high, low = value
            if low <= wavenumber <= high:
                matched_groups.append(group)
        else: # Multiple ranges
            for high, low in value:
                if low <= wavenumber <= high:
                    matched_groups.append(group)

```

Functional group		Mode of vibration	Wavenumber (cm ⁻¹)	Intensity
C–H	Alkanes	(stretching)	3000-2850	Strong
	–CH ₃	(bending)	1450 and 1375	Medium
	–CH ₂ –	(bending)	1465	Medium
Alkenes	(stretching)	3100-3000	Medium	
	(out of plane bending)	1000-650	Strong	
Aromatics	(stretching)	3150-3050	Strong	
	(out of plane bending)	900-690	Strong	
Alkyne	(stretching)	ca. 3300	Strong	
	Aldehyde	2900-2800 2800-2700	Weak Weak	
C=C	Alkane	not interpretative useful		
	Alkene	1680-1600	Medium to Weak	
	Aromatic	1600 and 1475	Medium to Weak	
C≡C	Alkyne	2250-2100	Medium to Weak	

```

# Dictionary for functional groups
functional_groups = {
    "C–H Alkanes": (3000, 2850),
    "C–H Alkanes –CH3": [(1450, 1450), (1375, 1375)],
    "C–H Alkanes –CH2–": (1465, 1465),
    "C–H Alkenes": (3100, 3000),
    "C–H Aromatics": [(3150, 3050), (900, 690)],
    "C–H Alkyne": (3300, 3300),
    "C–H Aldehyde": [(2900, 2800), (2800, 2700)],
    "C=C Alkene": (1680, 1600),
    "C=C Aromatic": [(1600, 1600), (1475, 1475)],
    "C≡C Alkyne": (2250, 2100),
    "C=O Aldehyde": (1740, 1720),
    "C=O Ketone": (1725, 1705),
    "C=O Carboxylic acid": [(1725, 1700), (3400, 2400)],
    "C=O Ester": (1750, 1730),
    "C=O Amide": (1670, 1640),
    "C=O Anhydride": [(1810, 1810), (1760, 1760)],
    "C=O Acid chloride": (1800, 1800),
    "C–O Alcohols, Ethers, Esters": (1300, 1000),
    "O–H Alcohols, Phenols (Free)": (3650, 3600),
    "O–H Alcohols, Phenols (H-Bonded)": (3500, 3200),
    "N–H Primary and Secondary Amines": [(3500, 3100), (1640, 1550)],
    "C–N Amines": (1350, 1000),
    "C=N Imines and Oximes": (1690, 1640),
    "C≡N Nitriles": (2260, 2240),
    "X=C=Y Allenes, Ketenes, Isocyanates": (2270, 1950),
    "N=O Nitro (R–NO2)": [(1550, 1550), (1350, 1350)],
    "S–H Mercaptans": (2550, 2550),
    "S=O Sulfoxides": (1050, 1050),
    "C–X Fluoride": (1400, 1000),
    "C–X Chloride": (800, 600),
    "C–X Bromide and Iodide": (667, 0),
}

```

Adding the Values to the Table

```
# Style the highest peak row
if index == highest_peak_index:
    wavenumber = f"\033[1m\033[31m{wavenumber:.2f}\033[0m" # Bold + Red
    category = f"\033[1m\033[31m{category}\033[0m"           # Bold + Red

# Add row to the table
functional_groups_str = ', '.join(matched_groups) if matched_groups else "No matching functional group"
table_data.append([wavenumber, functional_groups_str, category])

# Print the combined table
print(tabulate(table_data, headers="firstrow", tablefmt="grid"))
print("\n")
```

Wavenumber	Functional Groups	Absorbance
2752.38	C-H Aldehyde, C=O Carboxylic acid	Medium
2768.93	C-H Aldehyde, C=O Carboxylic acid	Medium
2800.38	C-H Aldehyde, C=O Carboxylic acid	High
2822.57	C-H Aldehyde, C=O Carboxylic acid	High
2867.3	C-H Alkanes, C-H Aldehyde, C=O Carboxylic acid	Medium
2898.76	C-H Alkanes, C-H Aldehyde, C=O Carboxylic acid	Medium
2928.58	C-H Alkanes, C=O Carboxylic acid	Low
2954.4	C-H Alkanes, C=O Carboxylic acid	Low

Printed Results

Let users input the range of the wavenumber

Default range: Start Wavenumber = 587.63, End Wavenumber = 3902.94
If you wish to analyze a specific range, you can enter the range below:

Start Wave... 2500

End Waven... 3000


```
# Interactive analysis for a user-defined range
def interactive_analysis():
    interact(
        analyze_range,
        start=widgets.FloatText(value=data['Wavenumber'].min(), description='Start Wavenumber:'),
        end=widgets.FloatText(value=data['Wavenumber'].max(), description='End Wavenumber:')
    )


# Function to analyze the selected range
def analyze_range(start, end):
    combined_analysis(start, end)

# Option to analyze a specific range
print(f"Default range: Start Wavenumber = {data['Wavenumber'].min():.2f}, End Wavenumber = {data['Wavenumber'].max():.2f}")
print("If you wish to analyze a specific range, you can enter the range below:")
interactive_analysis()
```


RESULT

Formaldehyde







Wavenumber	Functional Groups	Absorbance
1505.08	No matching functional group	Low
1717.38	C=O Ketone, C=O Carboxylic acid	Medium
1748.83	C=O Ester	Medium
1771.02	No matching functional group	Medium
2752.38	C-H Aldehyde, C=O Carboxylic acid	Medium
2768.93	C-H Aldehyde, C=O Carboxylic acid	Medium
2800.38	C-H Aldehyde, C=O Carboxylic acid	High
2822.57	C-H Aldehyde, C=O Carboxylic acid	High
2867.3	C-H Alkanes, C-H Aldehyde, C=O Carboxylic acid	High
2898.76	C-H Alkanes, C-H Aldehyde, C=O Carboxylic acid	Medium
2928.58	C-H Alkanes, C=O Carboxylic acid	Low
2954.4	C-H Alkanes, C=O Carboxylic acid	Low




Phenol







Wavenumber	Functional Groups	Absorbance
508.768	C-X Bromide and Iodide	Low
689.234	C-X Chloride	Medium
749.41	C-H Aromatics, C-X Chloride	High
811.522	C-H Aromatics	Medium
1177.45	C-O Alcohols, Ethers, Esters, C-N Amines, C-X Fluoride	Medium
1217.91	C-O Alcohols, Ethers, Esters, C-N Amines, C-X Fluoride	Medium
1252.43	C-O Alcohols, Ethers, Esters, C-N Amines, C-X Fluoride	Medium
1344.13	C-N Amines, C-X Fluoride	Low
1475.34	No matching functional group	Medium
1500.93	No matching functional group	High
1598.63	N-H Primary and Secondary Amines	High
3355.18	C=O Carboxylic acid, O-H Alcohols, Phenols (H-Bonded), N-H Primary and Secondary Amines	Medium
3613.62	O-H Alcohols, Phenols (Free)	Low




Salicylic Acid





Wavenumber Functional Groups	Absorbance
3277 C=O Carboxylic acid, O-H Alcohols, Phenols (H-Bonded), N-H Primary and Secondary Amines	Low
3042 C-H Alkenes, C=O Carboxylic acid	Low
2881 C-H Alkanes, C-H Aldehyde, C=O Carboxylic acid	Low
1679 C=C Alkene, C=N Imines and Oximes	High
1625 C=C Alkene, N-H Primary and Secondary Amines	Medium
1496 No matching functional group	Medium
1479 No matching functional group	Medium
1458 No matching functional group	High
1302 C-N Amines, C-X Fluoride	Medium
1258 C-O Alcohols, Ethers, Esters, C-N Amines, C-X Fluoride	High
1217 C-O Alcohols, Ethers, Esters, C-N Amines, C-X Fluoride	Medium
1199 C-O Alcohols, Ethers, Esters, C-N Amines, C-X Fluoride	Medium
1160 C-O Alcohols, Ethers, Esters, C-N Amines, C-X Fluoride	Medium
769 C-H Aromatics, C-X Chloride	Medium
702 C-H Aromatics, C-X Chloride	Medium



Group Work

Charoline 黃菊花
B10611026

Made the plot and
finding the peaks

Asia 黃夏彤
B11502120

Adding peak description
and user interface

Teresa 陳翠珍
B11504086

Converting JDX file to
excel file

Thalia 郭伶美
B11B02043

Code fixing, evaluating,
and result testing



THANK YOU!